1. Intended use / Indication

Light-curable polymerizable resin intended to be used in conjunction with curing light equipment. Indicated for the fabrication of tough-elastic earmolds.

2. Contraindication

audioprint® GR-11 is contraindicated for every application that is not part of the indication (see above).

3. Patient Target Group

Persons for whom an earmold is to be created.

4. Intended Users

Earmold laboratory, ENT specialist, hearing aid acoustician

5. Requirements

Software:

Cyflex AG · Binzmühlstrasse 15, 8050 Zürich, Switzerland
3SHAPE A/S · Holmens Kanal 7, 1060 Copenhagen, Denmark

Hardware (3D Printing):

ASIGA - 2/19-21 Bourke Road, Alexandria, NSW 2015, Australia

Hardware (Post Curing):

pro3dure medical (CD-1, CD-2) · Am Burgberg 13, D-58642 Iserlohn, Germany

Additional specifications on www.pro3dure.com.

6. Material

Contains: Alkoxilated bisphenole-A-dimethacrylate, initiators, dyes, stabilisers, and pigments.

7. Geometric presetting

According to the appropriate design specifications (templates).

8. Material parameters

Radiation penetration depth controllable by exposure time.

9. Manufacturing process (fig. 1-10)

1. Prepare data (CAD & build preparation).
2. Choose process parameters (build style, etc.).
3. Transfer prepared data to 3D printing device.
4. Prepare 3D printing – shake bottle.
5. Fill resin tank of 3D printing device.
6. Build the parts.
7. Clean parts (with IPA ≥ 97 % or equivalent cleaning agent) approx. 5 min. in an ultrasonic bath or equivalent device – precleaning recommended.
8. Dry parts, until there are no residues of IPA or equivalent cleaning agent.
9. Post curing (5 min.): inert atmosphere recommended (use adequate light curing device).
10. Finish parts.

10. Finishing processes

Polishing

11. Notice


Please follow the instructions provided by the software manufacturer in terms of parameter settings and design recommendations. Please follow the instructions provided by the hardware manufacturer in terms of parameter settings/printing & post curing recommendations. To avoid detrimental effects on material quality do not expose the liquid material to irradiation under any circumstances. Deviations from the described manufacturing processes or storing conditions may lead to different mechanical and optical properties of the material. Ensure personal protective gear during processing. According to the EU Medical Devices Regulation, users/patients are obliged to report serious events with a medical device to the manufacturer and to the competent authority of the country in which they occurred. **Caution:** Avoid breathing dust/fume/gas/mist/vapours/spray. Wash hands thoroughly after handling. Use only outdoors or in a well-ventilated area. Contaminated work clothing should not be allowed out of the workplace. Wear protective gloves/protective clothing/eye protection/face protection. IF ON SKIN: Wash with plenty of soap and water. IF INHALED: Remove person to fresh air and keep comfortable for breathing. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. Call a POISON CENTER or doctor/physician if you feel unwell. Specific treatment (see information on this label). If skin irritation occurs: Get medical advice/attention. If skin irritation or rash occurs: Get medical advice/attention. If eye irritation persists: Get medical advice/attention. Take off contaminated clothing and wash it before reuse. Store in a well-ventilated place. Keep container tightly closed. Store locked up. Dispose of contents/container according to official regulations.

12. Hazard statements

Causes skin irritation. May cause an allergic skin reaction. Causes serious eye irritation. May cause respiratory irritation.

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Distributor (US):
pro3dure medical LLC
9825 Valley View Road · Eden Prairie, MN 55344
Phone: 952-426-1928 · Fax: 952-681-7515
info@pro3dure.com · www.pro3dure.com · Made in Germany

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Hardware (3D Printing):
ASIGA · 2/19-21 Bourke Road, Alexandria, NSW 2015, Australien

Hardware (Post Curing):
pro3dure medical (CD-1, CD-2) · Am Burgberg 13, D-58642 Iserlohn, Deutschland

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