1. Indication
printodont[®] GR-14.1 denture from pro3dure is a light-curable poly-merizable resin intended to be used in conjunction with extraoral curing light equipment. **printodont[®] GR-14.1 denture** is indicated for the fabrication and repair, by additive manufacturing, of full and partial removable dentures and baseplates.

2. Contraindication
printodont[®] GR-14.1 denture restorations are contraindicated ...

1. ... if a patient is known to be allergic to any of the ingredients.
2. ... direct intraoral contact of uncured or partially uncured material.
3. ... for every application that is not part of the indication (see above).

3. Requirements

File Format
 STL (see geometric presettings for CAD)

Hardware (3D Printing):

ASIGA · 2/19-21 Bourke Road, Alexandria, NSW 2015, Australia
 Structo Pte Ltd · 60 Mac Pherson Road, #06-08A, The Siemens Centre, Singapore (348615)
 MIICRAFT · 27 Queen St. East Suite 1401 Toronto, Ontario M5C 2M6, Canada
 DAZZ3D · Rm K-2F, Building A, No. 9 Zhongshanyuan Rd., Shenzhen, China 518052

Hardware (Post Curing):

pro3dure medical (CD-1, CD-2) · Am Burgberg 13, D-58642 Iserlohn, Germany
 Structo Pte Ltd · 60 Mac Pherson Road, #06-08A, The Siemens Centre, Singapore (348615)
 NK Optik (Ottoflash G171) · Isarstraße 2, 82065 Baierbrunn, Germany
 Prusa Research a.s. (Prusa CW1) · Partyzánská 188/7a, 170 00 Praha 7, Czech Republic

Bonding agent for fixing denture teeth in the denture base:

VITA Zahnfabrik H. Rauter GmbH & Co. KG · Spitalgasse 3, D-79713 Bad Säckingen, Germany

4. Material

printodont[®] GR-14.1 denture consists of functional (meth)acrylic resins, initiators, dyes and stabilisers.

5. Geometric presetting

Lingual Ridge ≥ 2.5 mm (Maxilla / Mandibula)
Palatal/Lingual ≥ 3.0 mm (Maxilla) / ≥ 2.5 (Mandibula)
Facial/Buccal ≥ 2.5 mm (Maxilla / Mandibula)

6. Material parameters

Radiation penetration depth controllable by exposure time
 ≥ 50 µm (4 mils)
 ≥ 100 µm (5 mils)

7. Manufacturing process (fig. 1-10, page 1)

1. Prepare data (CAD & build preparation).
2. Choose process parameters (build style, etc.).
3. Transfer prepared data to 3D printing device.
4. Prepare 3D printing – shake bottle.
5. Fill resin tank of 3D printing device.
6. Build the parts.
7. Clean parts (with IPA ≥ 97 % or equivalent cleaning agent) approx. 4 min. in an ultrasonic bath or equivalent device – precleaning recommended.
8. Dry parts, until there are no residues of IPA or equivalent cleaning agent.
9. Post curing (10 min.): inert atmosphere recommended (use adequate light curing device).
10. Finish parts.

8. Bonding of denture teeth

If using printed denture teeth – see additional IFU **printodont[®] GR-17.1 temporary It** for details. In order to realize the optimal bonding of artificial teeth to the dental object, it can be necessary to roughen the surface of the denture and to use an appropriate adhesive. The suitability of a bonding agent or adhesive must be approved prior use. Recommended bonding agent is VITA VIONIC[®] BOND. The commonly used techniques and dental work steps must be carefully considered when manufacturing the dentures. Please follow the instructions and indications of the corresponding suppliers.

9. Notice

Please follow the instructions provided by the software manufacturer in terms of parameter settings and design recommendations. Please follow the instructions provided by the hardware manufacturer in terms of parameter settings/printing & post curing recommendations. To avoid detrimental effects on material quality do not expose the liquid material to irradiation under any circumstances. Deviations from the described manufacturing processes or storing conditions may lead to different mechanical and optical properties of the material. Ensure personal protective gear during processing.


Caution: Polymerized resins are chemically resistant – avoid stains on clothing! Avoid any contact with skin and eyes! In case of accidental contact, rinse with adequate running water, consult a doctor if necessary. The lot number and the best before date are indicated on each packaging. In case of claims please always indicate the lot number of the product. Do not use the product after expiry of the best before date. Dispose of contents/container according to official regulations.

10. Hazard statements

Causes skin irritation. May cause an allergic skin reaction. Causes serious eye irritation. May cause respiratory irritation.

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1 kg:
 λ ≤ 405 nm

orange-pink/
 orange-rose/
 rosa-naranja/
 rosa-arancione
 REF: D1001302

light-pink/
 rose-clair/
 rosa-chiaro/
 rosa-chiaro
 REF: D1001303

deep-pink/
 dunkel-pink/
 rose-profond/
 rosa-profundo/
 rosa-scuro
 REF: D1001304

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Distributor (US):
 pro3dure medical LLC
 9825 Valley View Road · Eden Prairie, MN 55344
 Phone: 952-426-1928 · Fax: 952-681-7515
 info@pro3dure.com · www.pro3dure.com · Made in Germany

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