1. Intended use / Indication
otosil® IS-1.1 is an addition-curing 2K silicone for molding the external auditory canal (ear impression).

2. Contraindication

otosil® IS-1.1 ear impression silicone is contraindicated ...
1. ... if it is known that a patient is allergic to one of the ingredients.
2. ... for every application that is not part of the indication (see above).

3. Patient Target Group

Persons for whom an earmold is to be created.

4. Intended Users

Earmold laboratory, ENT specialist, hearing aid acoustician

5. Requirements

Impressions of the ear canal may only be carried out by trained specialists (hearing aid acousticians or ENT specialists). The work instructions and precautionary measures must be strictly observed. Failure to do so can result in irreparable damage to the ear or eardrum.

6. Material

otosil® IS-1.1 is an addition-curing 2K silicone material with special flow properties for pressure-free precision molding of the external auditory canal.

7. Preparations, general tips

- Careful preparation as well as thorough cleaning of the ear canal before taking the impression ensure better and more reliable results.
- Ensure that all traces of impression material are removed from the patient's ear.
- The hardened impression material otosil® IS-1.1 is chemically inert. Avoid staining clothes.
- We recommend standard nitrile or polyethylene gloves.

8. Prepare the ear

The ear canal and eardrum must be examined before the impression is taken. Ear impressions should not be taken if there are indications of characteristic features (e.g. perforation or inflammation of the tube or eardrum). The ear canal needs to be cleaned and depilated. In the event of insufficient cleaning, the hardening mechanism of the impression material can be impaired (e.g. by cerumen). In addition to cleaning, it is recommended to create a thin separating layer between the impression material and the entire skin contact area (e.g. ear canal, concha and edge areas of the auricle). To do this, the relevant areas must be moistened with a lubricating solution using a cotton ball or cotton pad. By conditioning the ear in the aforementioned way, hardened impression material is more easily removed and a more homogeneous hardened surface of the ear impression (reduced inhibition layer) is obtained.

Before the impression is taken, the eardrum must always be protected with an impression plug at the end of the external ear canal. The impression silicone otosil® IS-1.1 is now applied directly into the conditioned ear of the patient as described under point 9.

9. Mixing and dosing (fig. 1-5)

- Before the impression is taken, the eardrum must always be protected with an impression plug at the end of the external ear canal (fig. 1).
- Insert the cartridge (fig. 2) into the dispensing gun.
- Unscrew the cap. To ensure proper flow and mixing from both ports, carefully extrude a small amount of the material (fig. 3).
- Insert the static mixer into the guided grooves of the cartridge. Turn the mixing cannula counterclockwise as far as it will go (fig. 4). If necessary, mount an intra-tip. As an alternative, static mixers with a larger diameter (9 mm) can be used to reduce the forces during dispensing. The dispenser is now ready to mix the 2-component impression silicone in any required amount.
- Inject the silicone material slowly into the ear canal until the concha is completely filled (fig. 5). The special flow properties of the imprint silicone otosil® IS-1.1 from pro3dure enable an almost pressure-free, constant flow into the ear. Always keep the end of the static mixer on the surface of the injected material during injection. To seal the cartridge, leave the static mixer on after use. After curing, carefully remove the ear impression and check the ear again.

10. Follow-up treatment and control

Every ear impression carried out must be carefully completed by trained specialists as part of a follow-up examination and control. In particular, the auditory canal must be assessed.

11. Note

The impression of the external auditory canal may only be carried out by trained specialists. The work instructions and precautionary measures must be strictly observed. Failure to do so can result in irreparable damage to the ear or eardrum. pro3dure medical is not liable for damage caused by improper use of the impression material.

According to the EU Medical Devices Regulation, users/patients are obliged to report serious events with a medical device to the manufacturer and to the competent authority of the country in which they occurred.


Despite the fact that the high biocompatibility of silicone materials has been proven over decades in exceptional cases, undesirable reactions of the immune system such as allergies cannot be absolutely ruled out. If in doubt, we recommend that you contact your doctor before using the material and carry out an allergy test.

12. Hazard warnings

May be harmful to aquatic organisms with long-term effects. Disposal in accordance with official regulations.

Impression silicone IS-1.1

Instructions for use

Physical properties*:

otosil® IS-1.1

- Shore hardness A (ISO 48-4): 45 ± 5
- Mixing volume: 50 ml Cartridge
- Mixing ratio: 1:1
- Setting time: < 3 min. at body temperature
- Storage: 18-28 °C

Ordering information:

otosil® IS-1.1

Standard packing:

Cartridge
8 x 50 ml
(à 2 x 25 ml A + B)

Base:
orange-opaque

+

Catalyst:
white

REF: A0010100

Cannula:

S System
4,2 mm
yellow
50 Pieces

REF: A4020141

Dispenser:
DS50

REF: A4020139

* These data come from measurements of a representative sample that was determined as part of our quality assurance.

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