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
SysTherm® PU is a 2-component polymer for the production of earmolds with thermoeffect. The cured **SysTherm® PU** material is characterized by an initial Shore A hardness of 85 at room temperature (23 °C) which changes to 65 Shore A at body temperature (37 °C). Accordingly, **SysTherm® PU** earmolds can be worked on at room temperature. Moreover, the end customer can precisely position a **SysTherm®** earmold in their ear canal. At body temperature the earpiece becomes soft and flexible, generating outstanding wearing comfort in combination with a perfect sealing during jaw movements. **SysTherm®**s excellent recovery ability ensures that the earmold regains its original shape after deformation. When using a digital workflow the **SysTherm® PU** material can be processed with 3D printed castforms (pro3dure's **audioprint® GR 12.2**) and the **SysTherm® SPT-1** Separator or solely with pro3dure's **SysTherm® Cast**. Additionally, the usage of the material in the traditional PnP process with a silicone negative form (pro3dure's **FS-2**) is also possible. Due to its chemical base **SysTherm®** is moisture sensitive and should be stored dry at room temperature (23 °C).

- 1) Open the **SysTherm®** secondary packaging (fig. 1).
- 2) Twist off the cartridge top and attach the mixing canule (fig. 2).
- 3) Insert the cartridge into the **SysTherm®** dispenser (fig. 3).
- 4) Before dispensing into a mold please pump a pearl-sized amount of **SysTherm® PU** through the mixing tip in order to ensure a perfect mixture of the material (fig. 4).
- 5) Depending on the chosen production process please proceed to one of the corresponding sections below.
- 6) For storing of an opened **SysTherm® PU** cartridge please leave the mixing canula on the cartridge (fig. 6).

A)


Production process based on 3D-printed molds (**audioprint® GR 12.2**) with cast separator (**SysTherm® SPT-1** Separator):

- 1) Printing cast shells out of **audioprint® GR 12.2** (Make sure air-vents are the same height as filling inlet, otherwise PU material will leak out during curing) (fig. 7).
- 2) After printing thoroughly clean the parts using IPA and then fresh IPA, for 5 min. each (fig. 8).
- 3) Dry the parts with compressed air (fig. 9).
- 4) Dry the parts in oven at 40 °C for at least 30 min. to eliminate residual IPA. Parts must be odor free.
- 5) Cure the parts with nitrogen for 2 min. in **CD-2**.
- 6) Dip the casts into separator liquid **SysTherm® SPT-1** (fig. 10).
- 7) Thoroughly dry the cast with compressed air (fig. 11).
- 8) Dry in oven at 60 °C for 120 min. (fig. 14).
- 9) Thoroughly dry the cast with compressed air.
Note: the separator will stay as a thin liquid film on the surface (fig. 11).
- 10) Fill in **SysTherm® PU** carefully and avoid trapped air. Fill the inlet as a reservoir for possible polymerization shrinkage (fig. 13).
- 11) Put the filled casts into a **dry pressure pot** – make sure **SysTherm®** will not flow out.
- 12) Cure **SysTherm®** at 60 °C / 2 bar for 2 hours. Keep heating up time into account (fig. 14).
- 13) Cool down the filled parts and carefully crack the shells. Use eye protection and appropriate tools (fig. 15).
- 14) Peel off cracked cast parts.
- 15) Clean and finish **SysTherm®** parts.

Indications for use:
silicone earmould
material

Technical data:

- Mixing ratio: 4:1
- Setting time: 120 min. at 60 °C / 140 °F in pressure pot
- Final hardness: Shore A: 23:85 ±5
Shore A: 37:70 ±5
- Application: at 23 ± 2 °C, 50 ± 5 % rel. humidity
- Storage:

Ordering information:

SysTherm® PU
Polymer
2 x 50 ml
4 mixing canulas
clear-transparent
item no.: A4020124
rose-transparent
item no.: A4020144
reddish-orange
item no.: A4020143

SysTherm® Cast
1 kg
item no.: A4020125


SysTherm® SPT-1
500 g
item no.: A4020126

SysTherm®
mixing canulas
25 pieces
item no.: A4020127

SysTherm®
4:1 dispenser
1 piece
item no.: A4020128

CE
0044

These data result from measurements of a representative sample, which were determined within the scope of our quality assurance.

B)


Production Process based on 3D-printed molds (**SysTherm®**):

- 1) Printing cast shells out of **SysTherm® Cast**. Make sure the air-vents are the same height as filling inlet, otherwise PU material will leak out during curing (fig. 7).
- 2) After printing thoroughly dry the parts with compressed air to remove residual liquid resin (fig. 9).
- 3) Cure the parts with nitrogen for 2 min. in **CD-2**.
- 4) Do NOT use any kind of cast separators.
- 5) Dry the casts in oven at 40 °C for 60 min. to remove any moisture (fig. 12).
- 6) Fill in **SysTherm® PU** carefully and avoid trapped air. Fill the inlet as a reservoir for possible polymerization shrinkage (fig. 13).
- 7) Put the filled casts into a dry pressure pot. Make sure **SysTherm®** will not flow out.
- 8) Cure **SysTherm®** at 60 °C for 2 hours. Keep heating up time into account (fig. 14).
- 9) Place the casts into a beaker with hot water (60-70 °C) and wait 30 min. Cast material will swell (fig. 16).
- 10) Peel off swelled cast material and remove parts (fig. 15).
- 11) Clean and finish **SysTherm®** parts.

C)

Production process based on PNP method using form silicones :

- 1) Preparation of negative mould out of **FS-2** silicone (see instruction for use).
- 2) If necessary cleaning of mould and drying (60 °C, 120 min.) to avoid moisture (fig. 12).
- 3) Fill in **SysTherm® PU** carefully and avoid trapped air (fig. 18).
- 4) Put the filled casts into a dry pressure pot – make sure **SysTherm®** will not flow out.
- 5) Cure **SysTherm®** at 60 °C for 2 hours. Keep heating up time into account (fig. 14).
- 6) Clean and finish **SysTherm®** parts.

Temperature
limit

CE mark

Consult instructions
for use

Use by
date

Batch code

Catalogue
number

Date of
manufacture

idea to product.

Made in Germany
pro3dure medical GmbH
Am Burgberg 13
D - 58642 Iserlohn
Tel. +49 (0)2374.920050-0
Fax +49 (0)2374.920050-50
info@pro3dure.com
www.pro3dure.com