1. Intended use / Indication

The radiation-curing signature gel **SGNR-1** can be used advantageously to further individualize your earmolds. Due to its optimized thixotropy, the finest details of e.g. laser engravings can be filled without the material smearing or running during the working process.

2. Contraindication

SGNR-1 is contraindicated ...

- ... if a patient is known to be allergic to any of the ingredients.
- ... for every application that is not part of the indication (see above).

3. Patient Target Group

Individuals receiving hearing aid treatment.

4. Intended Users

ENT physician, hearing aid acoustician

5. Requirements

Hardware (Post Curing):

pro3dure medical (CD-1, CD-2) · Am Burgberg 13, D-58642 Iserlohn, Germany

6. Material

SGNR-1 consists of functional (meth)acrylic resins, initiators, dyes and stabilisers.

7. Geometric presetting

Not applicable.

8. Material parameters

Not applicable.

9. Manufacturing process (fig. 1-5)

1. Degrease the areas to be roughened with a suitable agent (e.g. isopropyl alcohol).
2. Apply **SGNR-1** gel to the engravings.
3. Carefully wipe away excess material.
4. Curing of the **SGNR-1** gel with suitable post-curing equipment (2 min. **CD-2**, under nitrogen).
5. Post-processing and finishing of the prepared area.

10. Post-processing and finishing of the repair area

When finishing the prepared areas, use the usual techniques and work steps (sanding/vibratory grinding). Please observe the instructions and notes of the relevant suppliers. Surface sealing can be done by polishing or painting (**L-1**).

11. Notice

Please follow the instructions provided by the hardware manufacturer in terms of parameter settings / post curing recommendations. To avoid detrimental effects on material quality do not expose the liquid material to irradiation under any circumstances. Deviations from the described manufacturing processes or storing conditions may lead to different mechanical and optical properties of the material. Ensure personal protective gear during processing.


Caution: Obtain special instructions before use. Do not handle until all safety precautions have been read and understood. Wear protective gloves / protective clothing / eye protection / face protection. If exposed: Call a POISON CENTER or doctor/physician. Store locked up. Dispose of contents/container in accordance with official regulation.

12. Hazard statements

May damage the unborn child. Suspected of damaging fertility.

**Signature gel · Signaturgel
Gel signature · Gel de firma
Gel di firma
SGNR-1**

Instructions for use · Gebrauchsanweisung
Mode d'emploi · Instrucciones de uso
Istruzioni per l'uso

Physical properties* /
Physikalische Eigenschaften* /
Propriétés physiques* /
Propiedades físicas* /
Proprietà fisiche*:

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