**1. Indication**  
printo<sup>®</sup> GR-18.2 IB is a resin for the generative production of dental indirect bonding trays based on image projection systems.

**2. Contraindication**  
printo<sup>®</sup> GR-18.2 IB is contraindicated for every application that is not part of the indication (see above).

**3. Requirements**  
Software – Information available from:  
exocad GmbH · Julius-Reiber-Straße 37, 64293 Darmstadt, Germany

**Hardware (3D Printing) – Information available from:**  
ASIGA · 2/19-21 Bourke Road, Alexandria, NSW 2015, Australia

**Hardware (Post Curing) – Information available from:**  
pro3dure medical (CD-2) · Am Burgberg 13, D-58642 Iserlohn, Germany

**4. Material**  
printo<sup>®</sup> GR-18.2 IB consists of functional (meth)acrylic resins.

**5. Geometric presetting**  
According to the software design specifications.

**6. Material parameters**  
Radiation penetration depth controllable by exposure time  
50 µm  
100 µm

**7. Manufacturing process (fig. 1–10, page 1)**

1. Prepare data (CAD & build preparation).
2. Choose process parameters (build style, etc.).
3. Transfer prepared data to 3D printing device.
4. Prepare 3D printing – shake bottle.
5. Fill resin tank of 3D printing device.
6. Build the parts.
7. Clean parts (with IPA ≥ 97 % or equivalent cleaning agent) approx. 4 min. in an ultrasonic bath or equivalent device – precleaning recommended.
8. Dry parts (with compressed air until there are no residues of IPA or equivalent cleaning agent).
9. Post curing (10 min.): inert atmosphere recommended (use adequate light curing device).
10. Finish parts.


**8. Finishing processes**  
Not applicable

**9. Notice**  
Please follow the instructions provided by the software manufacturer in terms of parameter settings and design recommendations. Please follow the instructions provided by the hardware manufacturer in terms of parameter settings/printing & post curing recommendations. To avoid detrimental effects on material quality do not expose the liquid material to irradiation under any circumstances. Deviations from the described manufacturing processes or storing conditions may lead to different mechanical and optical properties of the material. Ensure personal protective gear during processing.  
**Caution:** Polymerized resins are chemically resistant – avoid stains on clothing! Avoid any contact with skin and eyes! In case of accidental contact, rinse with adequate running water, consult a doctor if necessary. The lot number and the best before date are indicated on each packaging. In case of claims please always indicate the lot number of the product. Do not use the product after expiry of the best before date. Dispose of contents/container according to official regulations.

**10. Hazard statements**  
May cause an allergic skin reaction. Harmful if inhaled. Toxic to aquatic life with long lasting effects.

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**Distributor (US):**  
pro3dure medical LLC  
9825 Valley View Road · Eden Prairie, MN 55344  
Phone: 952-426-1928 · Fax: 952-681-7515  
info@pro3dure.com · www.pro3dure.com · Made in Germany

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