1. Indication
pro3dure's generative resin **printo[®] GR-10 guide** is a resin for the generative production of dental objects (drill-guides, splints) based on image projection systems.

2. Contraindication
printo[®] GR-10 guide is contraindicated ...

- ... if a patient is known to be allergic to any of the ingredients.
- ... for bruxism splints and splints with wall thickness ≤ 1.5mm.
- ... for denture bases.
- ... for every application that is not part of the indication (see above).

3. Requirements

Software – Information available from:

DENTAL WINGS INC - 160 Rue St-Viateur E Suite 710, Montréal QC H2T 1A8 CANADA
3SHAPE A/S - Holmens Kanal 7, 1060 Copenhagen, Denmark

Hardware (3D Printing) – Information available from:

ASIGA - 2/19-21 Bourke Road, Alexandria, NSW 2015, Australia
Structo Pte Ltd - 60 Mac Pherson Road, #06-08A, The Siemens Centre, Singapore (348615)
MIICRAFT - 27 Queen St. East Suite 1401 Toronto, Ontario M5C 2M6, Canada

Hardware (Post Curing) – Information available from:

pro3dure medical (CD-1, CD-2) - Am Burgberg 13, D-58642 Iserlohn, Germany
Structo Pte Ltd - 60 Mac Pherson Road, #06-08A, The Siemens Centre, Singapore (348615)
NK Optik (Otoflash G171) - Isarstraße 2, 82065 Baierbrunn, Germany

4. Material

printo[®] GR-10 guide consists of functional methacrylic resins.

5. Geometric presetting

Minimum wall thickness: 1.5 mm

6. Material parameters

Radiation penetration depth controllable by exposure time

- 25 µm
- 50 µm
- 100 µm

7. Manufacturing process (fig. 1–10, page 1)

- Prepare data (CAD & build preparation).
- Choose process parameters (build style, etc.).
- Transfer prepared data to 3D printing device.
- Prepare 3D printing – shake bottle.
- Fill resin tank of 3D printing device.
- Build the parts.
- Clean parts (with IPA ≥ 97 % or equivalent cleaning agent) approx. 4 min. in an ultrasonic bath or equivalent device – precleaning recommended.
- Dry parts (with compressed air until no residues of IPA or equivalent cleaning agent).
- Post curing (4 min.): inert atmosphere recommended (use adequate light curing device).
- Finish parts.

8. Finishing processes

Polishing
Coating (pro3dure L-1 recommended)

9. Notice

Please follow the instructions provided by the software manufacturer in terms of parameter settings and design recommendations. Please follow the instructions provided by the hardware manufacturer in terms of parameter settings/printing & post curing recommendations. To avoid detrimental effects on material quality do not expose the liquid material to irradiation under any circumstances. Deviations from the described manufacturing processes or storing conditions may lead to different mechanical and optical properties of the material. Ensure personal protective gear during processing.


Caution: Polymerized resins are chemically resistant – avoid stains on clothing! Avoid any contact with skin and eyes! In case of accidental contact, rinse with adequate running water, consult a doctor if necessary. The lot number and the best before date are indicated on each packaging. In case of claims please always indicate the lot number of the product. Do not use the product after expiry of the best before date. Dispose of contents/container according to official regulations.

10. Hazard statements

Causes skin irritation. May cause an allergic skin reaction. Causes serious eye irritation. May cause respiratory irritation.

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Structo Pte Ltd - 60 Mac Pherson Road, #06-08A, The Siemens Centre, Singapur (348615)