**1. Indication**  
printodent<sup>®</sup> GR-22 flex from pro3dure is intended for the production of orthodontic and dental objects (mouthguards, nightguards, snoring appliances and splints).

**2. Contraindication**  
printodent<sup>®</sup> GR-22 flex is contraindicated ...

- ... if a patient is known to be allergic to any of the ingredients.
- ... for denture bases.
- ... for every application that is not part of the indication (see above).

**3. Requirements**

**Software:**  
DENTAL WINGS INC · 160 Rue St-Viateur E Suite 710, Montréal QC H2T 1A8, Canada  
exocad GmbH · Julius-Reiber-Straße 37, 64293 Darmstadt, Germany  
3SHAPE A/S · Holmens Kanal 7, 1060 Copenhagen, Denmark

**Hardware (3D Printing):**

ASIGA · 2/19-21 Bourke Road, Alexandria, NSW 2015, Australia

**Hardware (Post Curing):**

pro3dure medical (CD-1, CD-2) · Am Burgberg 13, D-58642 Iserlohn, Germany

**4. Material**

printodent<sup>®</sup> GR-22 flex consists of functional (meth)acrylic resins.

**5. Geometric presetting**

Minimum wall thickness: 1.5 mm

**6. Material parameters**

Radiation penetration depth controllable by exposure time  
50 µm  
100 µm

**7. Manufacturing process (fig. 1-10, page 1)**

- Prepare data (CAD & build preparation).
- Choose process parameters (build style, etc.).
- Transfer prepared data to 3D printing device.
- Prepare 3D printing – shake bottle.
- Fill resin tank of 3D printing device.
- Build the parts.
- Clean parts (with IPA ≥ 97 % or equivalent cleaning agent) approx. 10 min. in an ultrasonic bath or equivalent device – precleaning recommended.
- Dry parts, until there are no residues of IPA or equivalent cleaning agent.
- Post curing (2 x 5 min. / wall thickness > 2.5 mm: 2 x 10 min.); inert atmosphere recommended (use adequate light curing device).
- Finish parts.

**8. Finishing processes**

Polishing

**9. Notice**

Please follow the instructions provided by the software manufacturer in terms of parameter settings and design recommendations. Please follow the instructions provided by the hardware manufacturer in terms of parameter settings/printing & post curing recommendations. To avoid detrimental effects on material quality do not expose the liquid material to irradiation under any circumstances. Deviations from the described manufacturing processes or storing conditions may lead to different mechanical and optical properties of the material. Ensure personal protective gear during processing.

According to the EU Medical Devices Regulation, users/patients are obliged to report serious events with a medical device to the manufacturer and to the competent authority of the country in which they occurred.


**Caution:** Polymerized resins are chemically resistant – avoid stains on clothing! Avoid any contact with skin and eyes! In case of accidental contact, rinse with adequate running water, consult a doctor if necessary. The lot number and the best before date are indicated on each packaging. In case of claims please always indicate the lot number of the product. Do not use the product after expiry of the best before date. Dispose of contents/container according to official regulations.

**10. Hazard statements**

Causes skin irritation. May cause an allergic skin reaction. Causes serious eye irritation. May cause respiratory irritation. Harmful to aquatic life with long lasting effects.

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light blue/  
hellblau/  
bleu clair/  
azul claro/  
azzurro  
REF: D1001703

500 g:  
λ ≤ 405 nm

clear-transparent/  
klar-transparent/  
incoloro-transparente/  
incoloro-trasparente/  
incoloro-transparente  
REF: D1001702

light blue/  
hellblau/  
bleu clair/  
azul claro/  
azzurro  
REF: D1001704

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