1. Intended use / Indication

GB The PRO3dent® splint I MSI blank is a disc intended for the fabrication of dental appliances such as mouthguards, nightguards, bruxism and TMJ splint appliances, and bite splints using CAD/CAM milling machines.

2. Contraindication

PRO3dent® splint I MSI is contraindicated ...

- 1. ... if a patient is known to be allergic to any of the ingredients.
- 2. ... for every application that is not part of the indication (see above).

3. Product Description

PRO3dent® splint I MSI includes a lactam-based, biomimetic solution for multi-species inhibition*. This reduces the formation of undesirable biofilms on orthodontic appliances without a microcidal effect.

*c. albicans / s. mutans

4. Patient Target Group

Persons being treated in the context of a dental procedure.

5. Intended Users

Dentist, dental technicians

6. Processing

- PRO3dent® splint I MSI milling blanks should be used only for the purpose of CAD/CAM production of dental occlusal splints. They can be used in all common CAD/CAM systems (fig. 1).
- Suitable cross-cut carbide cutters are to be used exclusively for the processing of the PRO3dent® splint I MSI material.
- PRO3dent® splint I MSI occlusal splints are designed and milled by CAD/CAM technology and qualified staff. After finishing, the thickness in the occlusal region of the splint must not be less than 0,9 mm. For aesthetic reasons, a labial reduction to 0,8 mm is permissible. Any corners and edges should be rounded off.
- PRO3dent® splint I MSI splints can be removed from the milling blank by using suitable cross-cut carbide cutters or cutting discs suitable for composites. Carefully cut through the retaining strips without pressure.
- In order to avoid plaque accumulation subsequent polishing of the splints is essential. Please pre-polish with suitable silicone polishers and goat's hair brushes. The high-luster polishing should be done with a corresponding polishing compound (fig. 2). As far as possible, please avoid heat generation during the polishing and finishing of the PRO3dent® splint I MSI splint. This guarantees the optimal fitting of the splint.
- **Please note:** Upon completion, use an ultrasonic bath with water at a maximum temperature of approx. 35 °C for cleaning. Cleaning solutions are not recommended.

* not applicable / nicht anwendbar / non applicable / non applicable / non applicable

** on the basis of the dimension of a specimen of 4 mm diameter and 1 mm height

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
PRO3dent® splint I MSI milling blank for the generation of occlusal splints

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9. Safety advice

- **Storage:** No special measures are required.
- **Disposal:** Dental objects made of PRO3dent® splint I MSI are not water soluble, are inert, and do not present a hazard for ground water. As such, they can be disposed of as domestic waste without any special precautions if not stated otherwise in your local disposal guidelines.
- **Warranty:** Our recommendations concerning the application technique, regardless of whether they are communicated in writing, orally, or by means of practical instructions, are based on our own experiences and tests. As such, they are intended as guidelines only. We are continually striving to improve our products. Consequently, we reserve the right to make changes to their design and composition.
- **Hazard Information:** Processing of PRO3dent® splint I MSI blanks produces dusts which can irritate the eyes, skin, and airways. As such, it is essential to ensure the protective gear at your workplace is in perfect working order.
- **pro3dure** is not liable for any damages caused by improper application of the material. To be used by trained specialist personnel for the purpose indicated only.

Physical properties ***/
Physikalische Eigenschaften ***/
Propriétés physiques ***/
Propiedades físicas ***/
Proprietà fisiche ***:

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